

# Recombinant Human Protein N-terminal Glutamine Amidohydrolase,NTAQ1

Catalog No: #AP76141

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Package Size: #AP76141-1 10ug #AP76141-2 50ug #AP76141-3 500ug #AP76141-4 1mg

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human Protein N-terminal Glutamine Amidohydrolase,NTAQ1
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q96HA8
Uniprot	Q96HA8
GeneID	55093;
Calculated MW	49.8kDa
Target Sequence	MSPILGYWKIKGLVQPTRLLEYLEEKYEEHLYERDEGDKWRNKKFELGLEFPNLPYYIDGDVKLTQSMAIIRYI ADKHNMLGGCPKERAIEISMLEGAVLDIRYGVSRAYSKDFETLKVDFLSKLPPEMLKMFEDRLCHKTYLNGDHV THPDFMLYDALDVVLYMDPMCLDAFPKLVCFKKRIEAIQIDKYLKSSKYIAWPLQGWQATFGGGDHPKSDL VPRGSMEGNGPAAVHYQPASPPRDACVYSSCYEENVWKLCEYIKNHDQYPLEECYAVFISNERKMPIWKQ QARPGDGPVIWDYHVLLHVSSGGQSFYDLDTVLPFPCLFDITYVEDAIKSDDDIHPQFRKFRVICADSYLKN FASDRSHMKDSSGNWREPPPPYPCJETGDSKMMLNDFISMDPKVGGAVYTLSEFTHRFGSKNC
Formulation	Supplied as a 0.2 µm filtered solution of PBS,100mM GSH,1% TritonX-100,15% glycerol,pH7.4.
Storage	Store at < -20C, stable for 6 months after receipt.Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only