

## Recombinant Mouse Interleukin-21,IL-21

Catalog No: #AP76142



Package Size: #AP76142-1 10ug #AP76142-2 50ug #AP76142-3 500ug #AP76142-4 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Mouse Interleukin-21,IL-21
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Mouse
Accession No.	Q9ES17
Uniprot	Q9ES17
GeneID	60505;
Calculated MW	14.4kDa
Target Sequence	MPDRLLIRLRHLIDIVEQLKIYENDLDPELLSAPQDVKGHCHEAAAFACFQKAKLKPSNPGNNKTFIIDLVAQLRR RLPARRGGKKQKHIACPCSDSYEKRTPKFLEFLERLKWLLQKMIHQHLS
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM PB,150mM NaCl,pH7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only