

Recombinant Mouse Interleukin-4, IL-4

Catalog No: #AP76145



Package Size: #AP76145-1 10ug #AP76145-2 50ug #AP76145-3 500ug #AP76145-4 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Mouse Interleukin-4, IL-4
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Mouse
Accession No.	P07750
Uniprot	P07750
GeneID	16189;
Calculated MW	13.4kDa
Target Sequence	MHGCDKNHLREIIGILNEVTGEGTPCTEMDVPNVLTATKNTTESELVCRASKVLRIFYLKHGKTPCLKKNSSVL MELQRLFRAFRCILDSSISCTMNESKSTSLKDFLESLSIMQMDYS
Formulation	Lyophilized from a 0.2 µm filtered solution of PBS, pH7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7C for 2-7 days. Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only