

Recombinant Human Coronin-6, CORO6 (N-6His)

Catalog No: #AP76147



Package Size: #AP76147-1 10ug #AP76147-2 50ug #AP76147-3 500ug #AP76147-4 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Coronin-6, CORO6 (N-6His)
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q6QEF8-4
Uniprot	Q6QEF8
GeneID	84940;
Calculated MW	28.3kDa
Target Sequence	MNHKVVHHHHHMPSPWGWFAVTACGAQRNNFEEPVALQEMDTSNGVLLPFYDPDSSIVYLCGKGDSSIRYF EITDEPPFVHYLNTFSSKEPQRGMGFMPKRGLDVSKEIARFYKLHERKCEPIIMTVPRKSDLFQDDLYPDTPG PEPALEADEWLSGQDAEPVLISLRDGYVPPKHRELRVTKRNILDVRPPSGPRRSQSASDAPLSQHTLETLEEI KALRERVQAQEQRITALENMLCELVDGTD
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM PB, 150mM NaCl, pH 7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7C for 2-7 days. Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only