

Recombinant Mouse IL-1 Receptor Antagonist Protein,IL-1ra,IL-1RA

Catalog No: #AP76149

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP76149-1 10ug #AP76149-2 50ug #AP76149-3 500ug #AP76149-4 1mg

Description

Product Name	Recombinant Mouse IL-1 Receptor Antagonist Protein,IL-1ra,IL-1RA
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Mouse
Accession No.	P25085
Uniprot	P25085
GeneID	16181;
Calculated MW	17.5kDa
Target Sequence	MRPSGKRPCMKQAFRIWDTNQKTFYLRNNQLIAGYLQGPNIKLEEKIDMVPIDLHSVFLGIHGGKLCCLSCAKSG DDIKLQLEEVNITDLSKNKEEDKRFTFIRSEKGPSTSFESAACPGWFLCTTLEADRPVSLTNTPEEPLIVTKFYF QEDQ
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM PB, 150mM NaCl, pH 7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only