

Recombinant Mouse  $\alpha$ -Synuclein, SNCA (N-6His)

Catalog No: #AP76152



Package Size: #AP76152-1 10ug #AP76152-2 50ug #AP76152-3 500ug #AP76152-4 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Mouse $\alpha$ -Synuclein, SNCA (N-6His)
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Mouse
Accession No.	O55042
Uniprot	O55042
GeneID	20617;
Calculated MW	15.9kDa
Target Sequence	MNHKVVHHHHHMDVFMKGLSKAKEGVVAAAETKQGVAAEAGKTKEGVLVYVGSKTKEGVVHGVTVAEKTKEQVTNVGGAVVTGVTAVAQKTVEGAGNIAAATGFVKKDQMGKGEEGYPQEGILEDMPVDPGSEAYEMPSEE GYQDYEPEA
Formulation	Lyophilized from a 0.2 $\mu$ m filtered solution of 20mM PB, 150mM NaCl, pH 7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7C for 2-7 days. Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only