

Recombinant Mouse Interleukin-2, IL2

Catalog No: #AP76153



Package Size: #AP76153-1 10ug #AP76153-2 50ug #AP76153-3 500ug #AP76153-4 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Mouse Interleukin-2, IL2
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Mouse
Accession No.	P04351
Uniprot	P04351
GeneID	16183;
Calculated MW	17.4kDa
Target Sequence	MAPTSSSTSSSTAEAAQQQQQQQQQQHLEQLLMDLQELLSRMENYRNLKLPRLMTFKFYLPKQATELKDL QCLEDELGPLRHVLDLTQSKSFQLEDAENFISNIRVTVVKLKGSDNTFECQFDDESATVVDFLRRWIAFCQSIIS TSPQ
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM Sodium Citrate, 5% Trehalose, pH4.5.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7C for 2-7 days. Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only