

# Recombinant Human Homeobox Protein Hox-B4,HOXB4,HOX-2F (N-6His)

Catalog No: #AP76154

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Package Size: #AP76154-1 10ug #AP76154-2 50ug #AP76154-3 500ug #AP76154-4 1mg

## Description

Product Name	Recombinant Human Homeobox Protein Hox-B4,HOXB4,HOX-2F (N-6His)
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P17483
Uniprot	P17483
GeneID	3214;
Calculated MW	29.8kDa
Target Sequence	MGSSHHHHHHSSGLVPRGSHMAMSSFLINSNYVDPKFPPEEYSQSDYLPSDHSPGYAGGQRRESSFQP EAGFGRRAACTVQRYAACRDPGPPPPPPPPPPPPPPGLSPRAPAPPAGALLPEPGQRCEAVSSSPPPPP CAQNPLHPSPSHSACEPVVYPWMRKVHVSTVNPNYAGGEPKRSRTAYTRQQVDKLLKEFHYNRYLTRRRR VEIAHALCLSERQIKIWFQNRMRKWKDKHKL PNTKIRSGGAAGSAGGPPGRPNGGPRAL
Formulation	Supplied as a 0.2 µm filtered solution of 4mM HCl.
Storage	Store at < -20C, stable for 6 months after receipt.Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only