

Recombinant Human Interleukin-17B,IL-17B (C-6His)

Catalog No: #AP76155



Package Size: #AP76155-1 10ug #AP76155-2 50ug #AP76155-3 500ug #AP76155-4 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Interleukin-17B,IL-17B (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q9UHF5
Uniprot	Q9UHF5
GeneID	27190;
Calculated MW	19.2kDa
Target Sequence	QPRSPKSKRKGQGRPGPLAPGPHQVPLDLVSRMKPYARMEEYERNIEEMVAQLRNSSELAQRKCEVNLQLW MSNKRSLSPWGYSHDPSRIPVDLPEARCLCLGCVNPFTMQEDRSMVSPVFSQVPVRRRLCPPPPRTGPC RQRAVMETIAVGCTCIFVDHHHHHH
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM PB, 150mM NaCl, pH 7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only