

# Recombinant Mouse Leukocyte Mono Ig-Like Receptor 1,LMIR1,CD300a (C-Fc)



Catalog No: #AP76156

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Package Size: #AP76156-1 10ug #AP76156-2 50ug #AP76156-3 500ug #AP76156-4 1mg

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Mouse Leukocyte Mono Ig-Like Receptor 1,LMIR1,CD300a (C-Fc)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Mouse
Accession No.	Q6SJQ0
Uniprot	Q6SJQ0
GeneID	217303;
Calculated MW	44.3kDa
Target Sequence	LHGPSTMSGSVGESLSVSCRYEEKFKTKDKYWCRVSLKILCKDIVKTSSEEARSGRVTIRDHPDNLFTVTYELSLTLEDADTYMCAVDISLFDGSLGFDKYFKIELSVVPSSEDPVSSPGPTLETTPVVSTSLPTKGPALGSNTEGHREHDYSQGLRVDDIEGRMDEPKSCDKTHTCPPCPAPELLGGPSVFLFPPKPKDTLMISRTPEVTCVWVDVSHEDPEVKFNWYVDGVEVHNAKTKPREEQYNSTYRVVSVLTVLHQDWLNGKEYKCKVSNKALPAPIEKTISKAKGQPREPQVYTLPPSREEMTKNQVSLTCLVKGFYPSDIAVEWESNGQPENNYKTTPVLDSDGSFFLYSKLTVDKSRWQQGNVFCFSVMHEALHNHYTQKSLSLSPGK
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM PB, 150mM NaCl, pH 7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only