

Recombinant Mouse IL-2R α ,IL-2RA,CD25 (C-Fc)

Catalog No: #AP76157



Package Size: #AP76157-1 10ug #AP76157-2 50ug #AP76157-3 500ug #AP76157-4 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Mouse IL-2R α ,IL-2RA,CD25 (C-Fc)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Mouse
Accession No.	P01590
Uniprot	P01590
GeneID	16184;
Calculated MW	55.5kDa
Target Sequence	ELCLYDPPEVPNATFKALSYKNGTILNCECKRGFRRLKELVYMRCLGNSWSSNCQCTSNSHDKSRKQVTAQL EHQKEQQTITDMQKPTQSMHQENLTGHCREPPPWKHEDSKRIYHFVEGQSVHYECIPGYKALQRGPAISICK MKCGKTGWTQPQLTCVDEREHRFLASEESQGSRNSSPESETSCPITTTDFPQPTETTAMTETFLVLTMEYKDI EGRMDEPKSCDKTHTCPPCPAPPELLGGPSVFLFPPKPKDTLMISRTPEVTCVVDVSHEDPEVKFNWYVDGV EVHNAKTKPREEQYNSTYRVVSVLTVLHQDWLNGKEYKCKVSNKALPAPIEKTISKAKGQPREPQVYTLPPSR EEMTKNQVSLTCLVKGFYPSDIAVEWESNGQPENNYKTPPVLDSDGSFFLYSKLTVDKSRWQQGNVFCSSV MHEALHNHYTQKSLSLSPGK
Formulation	Lyophilized from a 0.2 μ m filtered solution of PBS, pH7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only