

Recombinant Mouse Transforming Growth Factor β-1, TGFβ1, TGFB1

Catalog No: #AP76158

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP76158-1 10ug #AP76158-2 50ug #AP76158-3 500ug #AP76158-4 1mg

Description

Product Name	Recombinant Mouse Transforming Growth Factor β-1, TGFβ1, TGFB1
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Mouse
Accession No.	P04202
Uniprot	P04202
GeneID	21803;
Calculated MW	12.8kDa
Target Sequence	ALDTNYCFSSTEKNCCVRQLYIDFRKDLGWKWIHEPKGYHANFCLGPCPYIWSLDTQYSKVLALYNQHNP GASAPCCVPQALEPLPIVYYVGRKPKVEQLSNMIVRSCKCS
Formulation	Lyophilized from a 0.2 μm filtered solution of 50mM Glycine, 150mM NaCl, pH 2.5.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7C for 2-7 days. Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only