

Recombinant Mouse Serum Amyloid P Component, SAP, Pentraxin 2 (C-6His)



Catalog No: #AP76163

Orders: order@signalwayantibody.com

Package Size: #AP76163-1 10ug #AP76163-2 50ug #AP76163-3 500ug #AP76163-4 1mg

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Mouse Serum Amyloid P Component, SAP, Pentraxin 2 (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Mouse
Accession No.	P12246
Uniprot	P12246
GeneID	20219;
Calculated MW	24.9kDa
Target Sequence	QTDLKRKVFVFPRESETHVHKLIPHLEKPLQNFTLCFRTYSDLSRSQSLFSYSVKGRDNELLYKEKVGEYSLYI GQSKVTVRGMEEYLSPVHLCITWESSGIVEFWVNGKPVVKKSLQREYTVKAPPSIVLGQEQDNYGGGFQR SQSFVGEFSDLYMWDYVLTQDILFVYRDSPVNPINLNWQALNYEINGYVIRPRVWDVDHHHHHH
Formulation	Lyophilized from a 0.2 µm filtered solution of PBS, pH7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7C for 2-7 days. Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only