

Recombinant Mouse Activated Leukocyte Cell Adhesion Molecule,ALCAM,CD166 (C-Fc)



Catalog No: #AP76164

Orders: order@signalwayantibody.com

Package Size: #AP76164-1 10ug #AP76164-2 50ug #AP76164-3 500ug #AP76164-4 1mg

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Mouse Activated Leukocyte Cell Adhesion Molecule,ALCAM,CD166 (C-Fc)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Mouse
Accession No.	Q61490
Uniprot	Q61490
GeneID	11658;
Calculated MW	83.3kDa
Target Sequence	CRTSECCFQDPPYPDADSGSASGPRDLRCYRISSDRYECSWQYEGPTAGVSHFLRCLSSGRCCYFAAGSA TRLQFSDQAGVSVLYTVTLWVESWARNQTEKSPEVTLQLYNSVKYEPPLGDIKVSFLAGLRMEWETPDNQV GAEVQFRHRTPSSPWKLGDCGPQDDDTESCLCPLMNVAQEFQLRRRQLGSQGSSWSKWSSPVCVPPEN PPQPQVRFVSVEQLGQDGRRLTLKEQPTQLELPEGCQGLAPGTEVTYRLQLHMLSCPCAKATRTHLHGKM PYLSGAAYNVAVISSNQFGPGLNQTWHIPADTHTEPVALNISVGTNGTTMYWPARAQSMTYCIEWQPVGGQDG GLATCSLTAPQDPDPAGMATYSWSRESGAMGQEKCYYITIFASAHPEKLTWSTVLSTYHFGGNASAAGTPH HVSVKNHSLDSVSDWAPSLSTCPGVLKEYVVRCDREDSKQVSEHPVQPTETQVTLGSLRAGVAYTVQVR ADTAWLRGVWSQPQRFSIEVDDIEGRMDEPKSCDKTHTCPPCPAPELLGGPSVFLFPPKPKDTLMISRTPEVT CVVVDVSHEDPEVKFNWYVDGVEVHNAKTKPREEQYNSTYRVVSVLTVLHQDWLNGKEYKCKVSNKALPAPI EKTISKAKGQPREPQVYTLPPSREEMTKNQVSLTCLVKGFYPSDIAVEWESNGQPENNYKTTTPVLDSDGSFF LYSKLTVDKSRWQQGNVFCFSVMHEALHNHYTQKSLSLSPGK
Formulation	Lyophilized from a 0.2 µm filtered solution of PBS, pH7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only