

Recombinant Mouse Interferon γ Receptor 1,IFN- γ R1,CD119 (C-6His)

Catalog No: #AP76165

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP76165-1 10ug #AP76165-2 50ug #AP76165-3 500ug #AP76165-4 1mg

Description

Product Name	Recombinant Mouse Interferon γ Receptor 1,IFN- γ R1,CD119 (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Mouse
Accession No.	P15261
Uniprot	P15261
GeneID	15979;
Calculated MW	26.9kDa
Target Sequence	ALTSTEDPEPPSPVPPTNVLKSYNLPVVCWEYQNMSQTPIFTVQVKVYSGSWTDSCTNISDHCCNIYEQIMY PDVSAWARVKAKVGQKESDYARSKEFLMCLKGKVGPPGLEIRRKKEEQLSVLVHFPEVVVNGESQGTMTFGD GSTCYTFDYTVYVEHNRSGEILHTKHTVEKEECNETLCELNISVSTLDSRYCISVDGISSFQVRTEKSKDVCIP PFHDDRKDVDHHHHHH
Formulation	Lyophilized from a 0.2 μ m filtered solution of 20mM PB, 150mM NaCl, pH7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only