

Recombinant Human IL-12 Receptor Subunit β1,IL-12RB1,CD212 (C-Fc)

Catalog No: #AP76170

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP76170-1 10ug #AP76170-2 50ug #AP76170-3 500ug #AP76170-4 1mg

Description

Product Name	Recombinant Human IL-12 Receptor Subunit β1,IL-12RB1,CD212 (C-Fc)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P42701
Uniprot	P42701
GeneID	3594;
Calculated MW	84.2kDa
Target Sequence	<p>CRTSECCFQDPPYPDADSGSASGPRDLRCYRISSDRYECSWQYEGPTAGVSHFLRCLSSGRCCYFAAGSA TRLQFSDQAGVSVLYTVTLWVESWARNQTEKSPEVTLQLYNSVKYEPPLGDIKVSFLAGLRMEWETPDNQV GAEVQFRHRTSPSPWKLGDGCPQDDDTESCLCLEMNVAQEFQLRRRQLGSQGSSWSKWSSPVCVPPEN PPQPQVRFVSVEQLGQDGRRLTLKEQPTQLELPEGCQGLAPGTEVTYRLQLHMLSCPCAKATRTHLHGKM PYLSGAAYNVAVISSNQFGPGLNQTWHIPADTHTEPVALNISVGTNGTTMYWPARAQSMTYCIEWQPVGQDG GLATCSLTAPQDPDPAGMATYSWSRESGAMGQEKCYYITIFASAHPEKLTWSTVLSTYHFGGNASAAGTPH HVSVKNHSLDSVSDWAPSLSTCPGVLKEYVVRCDREDSKQVSEHPVQPTETQVTLGSLRAGVAYTVQVR ADTAWLRGVWSQPQRFSIEVDDIEGRMDEPKSCDKTHTCPPCPAPELLGGPSVFLFPPKPKDTLMISRTPEVT CVVVDVSHEDPEVKFNWYVDGVEVHNAKTKPREEQYNSTYRVVSVLTVLHQDWLNGKEYKCKVSNKALPAPI EKTISKAKGQPREPQVYTLPPSREEMTKNQVSLTCLVKGFYPSDIAVEWESNGQPENNYKTTTPVLDSDGSFF LYSKLTVDKSRWQQGNVFCFSVMHEALHNHYTQKSLSLSPGK</p>
Formulation	Lyophilized from a 0.2 μm filtered solution of PBS, pH7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7C for 2-7 days. Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only