

Recombinant Mouse OX40, TNFRSF4, CD134 (C-Fc)

Catalog No: #AP76173



Package Size: #AP76173-1 10ug #AP76173-2 50ug #AP76173-3 500ug #AP76173-4 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Mouse OX40, TNFRSF4, CD134 (C-Fc)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Mouse
Accession No.	P47741
Uniprot	P47741
GeneID	22163;
Calculated MW	48.1kDa
Target Sequence	VTARRLNCVKHTYPSGHKCCRECPGHGMVSRCDHTRDTLCHPCETGFYNEAVNYDTCKQCTQCNRHSGS ELKQNCPTQDQTVCRCPGTQPRQDSGYKLGVCVPCPPGHFSPGNNQACKPWTNCTLSGKQTRHPASDS LDAVCEDRSLATLLWETQRPTFRPTTVQSTTWPRTSELPSPPTLVTPGPDIEGRMDPKSCDKTHTCPPCP APELLGGPSVFLFPPKPKDTLMISRTPEVTCVVVDVSHEDPEVKFNWYVDGVEVHNAKTKPREEQYNSTYRV VSVLTVLHQDWLNGKEYKCKVSNKALPAPIEKTISKAKGQPREPQVYTLPPSRDELTKNQVSLTCLVKGFYPS DIAVEWESNGQPENNYKTPPVLDSDGSFFLYSKLTVDKSRWQQGNVFCSCVMHEALHNNHYTQKSLSLSPG K
Formulation	Lyophilized from a 0.2 µm filtered solution of PBS, pH7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7C for 2-7 days. Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only