

## Recombinant Human B7-1,CD80 (C-6His)

Catalog No: #AP76180



Package Size: #AP76180-1 10ug #AP76180-2 50ug #AP76180-3 500ug #AP76180-4 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human B7-1,CD80 (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P33681
Uniprot	P33681
GeneID	941;
Calculated MW	24.7kDa
Target Sequence	VIHVTKEVKEVATLSCGHNVSVEELAQTRIWQKEKMKVLTMMSGDMNIWPEYKNRTIFDITNNLSIVILALRPS DEGTYECVVLKYEKDAFKREHLAEVTLVKADFPPTPSISDFEIPTSNIRRIICSTSGGFPEPHLSWLENGEELNAI NTTVSQDPETELYAVSSKLDNFMTTNHSFMCLIKYGHRLRVNQTFNWNNTTKQEHPDNHHHHHHH
Formulation	Lyophilized from a 0.2 µm filtered solution of PBS, pH 7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7C for 2-7 days. Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only