

# Recombinant Human RANK L, TRANCE, TNFSF11 (N-6His)

Catalog No: #AP76182

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Package Size: #AP76182-1 10ug #AP76182-2 50ug #AP76182-3 500ug #AP76182-4 1mg

## Description

Product Name	Recombinant Human RANK L, TRANCE, TNFSF11 (N-6His)
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	O14788
Uniprot	O14788
GeneID	8600;
Calculated MW	22.4kDa
Target Sequence	MGSSHHHHHHSSGLVPRGSHMIRAEEKAMVDGSWLDLAKRSKLEAQPFAHLTINATDIPSGSHKVSLSWYHD RGWAKISNMTFSNGKLIVNQDGFYLYANICFRHHETSGDLATEYLQLMVYVTKTSIKIPSSHTLMKGGSTKYW SGNSEFHFYISINVGFFKLRSGEEISIEVSNPSLLDPDQDATYFGAFKVRDID
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM Tris, 150mM NaCl, pH8.0.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7C for 2-7 days. Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only