

Recombinant Human RANK, TNFRSF11A, CD265 (C-6His)

Catalog No: #AP76183

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP76183-1 10ug #AP76183-2 50ug #AP76183-3 500ug #AP76183-4 1mg

Description

Product Name	Recombinant Human RANK, TNFRSF11A, CD265 (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q9Y6Q6
Uniprot	Q9Y6Q6
GeneID	8792;
Calculated MW	21.1kDa
Target Sequence	IAPPCTSEKHYEHLGRCCNKCEPGKYMSSKCTTTSDSVCLPCGPDEYLD SWNEEDKCLLHKVCDTGKALVAV VAGNSTTPRRCACTAGYHWSQDCECCRRNTECAPGLGAQHPLQLNKD TVCKPCLAGYFSDAFSSTDKCRP WTNCTFLGKRVEHHGTEKSDAVCSSLPARKPPNEPHVYLPVDHHHHHH
Formulation	Lyophilized from a 0.2 µm filtered solution of PBS, pH 7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7C for 2-7 days. Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only