

Recombinant Mouse Neuroblastoma Suppressor of Tumorigenicity 1,NBL1 (C-6His)

Catalog No: #AP76185

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP76185-1 10ug #AP76185-2 50ug #AP76185-3 500ug #AP76185-4 1mg

Description

Product Name	Recombinant Mouse Neuroblastoma Suppressor of Tumorigenicity 1,NBL1 (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Mouse
Accession No.	Q61477
Uniprot	Q61477
GeneID	17965;
Calculated MW	18.4kDa
Target Sequence	<p> APPINKLALFPDKSAWCEAKNITQIVGHSGCEAKSIQNRACLGQCFSYSVPNTFFPQSTESLVHCDSCMPAQS MWEIVTLECPGHHEEVPRVDKLVKIVHCSCQACGKEPSHEGLNVYVQGEDSPGSQPGPHSHAHPPGGQTP EPEEPPGAPQVEEEGAEDVDHHHHHH </p>
Formulation	Lyophilized from a 0.2 µm filtered solution of PBS, pH 7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only