

Recombinant Mouse Tyrosine-Protein Kinase Receptor TYRO3,DTK (C-Fc)



Catalog No: #AP76187

Orders: order@signalwayantibody.com

Package Size: #AP76187-1 10ug #AP76187-2 50ug #AP76187-3 500ug #AP76187-4 1mg

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Mouse Tyrosine-Protein Kinase Receptor TYRO3,DTK (C-Fc)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Mouse
Accession No.	P55144
Uniprot	P55144
GeneID	22174;
Calculated MW	68.9kDa
Target Sequence	AGLKLMGAPVKMTVSQQQPVKLNCSVEGMEDPDIHWMKDGTVVQNASQVSISEHSWIGLLSLKSVERSDA GLYWCQVKDGEETKISQSVWLTVEGVPFFTVEPKDLAVPPNAPFQLSCEAVGPPEPVTIYWWRGLTKVGGPA PSPSVLNVGTQRTQTEFSCEARNIKGLATSRPAIVRLQAPPAAPFNNTVTTISSYNASVAWVPGADGLALLHSC TVQVAHAPGEWEALAVVVPVPPFTCLLRNLAPATNYSLRVRCANALGSPSYGDWVPFQTKGLAPARAPQNF HAIRTDGLILEWEEVIPEDPGEGPLGPYKLSWVQENGTQDELMVEGTRANLTDWDPQKDLILRVCASNAIGD GPWSQPLVSSHAGRQGGPPHSRTSVDDIEGRMDEPKSCDKTHTCPPCPAPELLGGPSVFLFPPKPKDTL MISRTPEVTCVVDVSHEDPEVKFNWYVDGVEVHNAKTKPREEQYNSTYRVVSVLTVLHQDWLNGKEYKCK VSNKALPAPIEKTISKAKGQPREPQVYTLPPSREEMTKNQVSLTCLVKGFYPSDIAVEWESNGQPENNYKTP PVLDSGDSFFLYSKLTVDKSRWQQGNVFCFSVMHEALHNHYTQKSLSLSPGK
Formulation	Lyophilized from a 0.2 µm filtered solution of PBS, pH 7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7C for 2-7 days. Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only