

# Recombinant Mouse SLAM Family Member 5,SLAMF5,CD84(C-6His)

Catalog No: #AP76191

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Package Size: #AP76191-1 10ug #AP76191-2 50ug #AP76191-3 500ug #AP76191-4 1mg

## Description

Product Name	Recombinant Mouse SLAM Family Member 5,SLAMF5,CD84(C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Mouse
Accession No.	Q18PI6
Uniprot	Q18PI6
GeneID	12523;
Calculated MW	23.8kDa
Target Sequence	KDADPVVMNGILGESVTFLNLIQEPKKIDNIAWTSQSSVAFIKPGVNKAEVTITQGTYKGRIEIIDQKYDLVIRDLR MEDAGTYKADINEENEETITKIYYLHIYRRLKTPKITQSLISSLNNTCNITLTCSEKEEKDVTYSWSPFGEKSNV LQIVHSPMDQKLYTCTAQNPVSNSSDSVTVQQPCTDTPSFHPRHAVLPVDHHHHHH
Formulation	Lyophilized from a 0.2 µm filtered solution of PBS, pH 7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only