

# Recombinant Mouse IGF Binding Protein 6,IGFBP-6 (C-6His)



Catalog No: #AP76192

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Package Size: #AP76192-1 10ug #AP76192-2 50ug #AP76192-3 500ug #AP76192-4 1mg

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Mouse IGF Binding Protein 6,IGFBP-6 (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Mouse
Accession No.	P47880
Uniprot	P47880
GeneID	16012;
Calculated MW	23.7kDa
Target Sequence	ALAGCPGCGAGMQTGCRGGCVVEEDAGSPADGCTEAGGCLRREGQPCGVYSPKCAPGLQCQPRENEEAP LRALLIGQGRCQRARGPSEETTKESKPQGGASRSRDTNHRDRQKNPRTSAAPIRPNPVQDSEMGPCRRHLD SVLQQLQTEVFRGGARGLYVPNCDLRGFYRKQQRSSQGNRRGPCWCVDPMGQPLVSPDGQGSTQCSA RSSGVDHHHHHH
Formulation	Lyophilized from a 0.2 µm filtered solution of PBS, EDTA,DTT,pH7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only