

Recombinant Mouse B7-2,CD86 (C-Fc)

Catalog No: #AP76198



Package Size: #AP76198-1 10ug #AP76198-2 50ug #AP76198-3 500ug #AP76198-4 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Mouse B7-2,CD86 (C-Fc)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Mouse
Accession No.	P42082
Uniprot	P42082
GeneID	12524;
Calculated MW	52.2kDa
Target Sequence	VETQAYFNGTAYLPCPFTKAQNISELVFWQDQQLVLYEHYLGTEKLDVNAKYLGRTSFDRNNWTLRLH NVQIKDMGSYDCFIQKKPPTGSILQQTLTELSVIANFSEPEIKLAQNVGTGNSGINLTCTSKQGHPKPKKMYFLIT NSTNEYGDNMQISQDNVTELFISISNLSLSPDPGVVHMTVVCVLETESMKISSKPLNFTQEFPSQTYWKEVD DIEGRMDEPKSCDKTHTCPPCPAPELLGGPSVFLFPPKPKDTLMISRTPEVTCVVVDVSHEDPEVKFNWYVD GVEVHNAKTKPREEQYNSTYRVVSVLTVLHQDWLNGKEYKCKVSNKALPAPIEKTISKAKGQPREPQVYTLPP SREEMTKNQVSLTCLVKGFYPSDIAVEWESNGQPENNYKTPPVLDSDGSFFLYSKLTVDKSRWQQGNVFS CSVMHEALHNHYTQKSLSLSPGK
Formulation	Lyophilized from a 0.2 µm filtered solution of PBS, pH 7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only