

# Recombinant Mouse S100 Calcium Binding Protein B,S100B (C-6His)



Catalog No: #AP76201

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Package Size: #AP76201-1 10ug #AP76201-2 50ug #AP76201-3 500ug #AP76201-4 1mg

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Mouse S100 Calcium Binding Protein B,S100B (C-6His)
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Mouse
Accession No.	P50114
Uniprot	P50114
GeneID	20203;
Calculated MW	11.8kDa
Target Sequence	MSELEKAMVALIDVFHQYSGREGDKHKLKKSSELKELINNELSHFLEEIKEQEVVDKVMETLDEDGDGECDFQE FMAFVAMVTTACHEFFEHELEHHHHHH
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM PB, 150mM NaCl, pH7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only