

Recombinant Mouse Interferon  $\alpha$ -2

Catalog No: #AP76202



Package Size: #AP76202-1 10ug #AP76202-2 50ug #AP76202-3 500ug #AP76202-4 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Mouse Interferon $\alpha$ -2
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Mouse
Accession No.	P01573
Uniprot	P01573
Calculated MW	19.5kDa
Target Sequence	MCDLPHTYNLRNKRALKVLAQMRRLPFLSCLKDRQDFGFPLEKVDNQQIQKAQAIPVLRDLTQQTLNLFSTKA SSAAWNATLLDSFCNDLHQQLNDLQTCLMQQVGVQEPPLTQEDALLAVRKYFHRITVYLREKKHSPCAWEVV RAEVWRALSSSVNLLPRLSEEKE
Formulation	Lyophilized from a 0.2 $\mu$ m filtered solution of 20mM PB, 100mM NaCl, pH7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7C for 2-7 days. Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only