

# Recombinant Mouse IL-1 Receptor-Like 1, IL-1RL1, IL-1 R4 (C-6His)

Catalog No: #AP76207

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Package Size: #AP76207-1 10ug #AP76207-2 50ug #AP76207-3 500ug #AP76207-4 1mg

## Description

Product Name	Recombinant Mouse IL-1 Receptor-Like 1, IL-1RL1, IL-1 R4 (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Mouse
Accession No.	P14719-2
Uniprot	P14719
GeneID	17082;
Calculated MW	36.6kDa
Target Sequence	SKSSWGLENEALIVRCPQRGRSTYPVEWYSDTNESIPTQKRNRIFVSRDRLKFLPARVEDSGIYACVIRSPNL NKTGYLNVTIHKKPPSCNIPDYLMYSTVRGSDKNFKITCPTIDLYNWTAPVQWFKNCKALQEPRFRAHRSYLF DNVTHDDEGDYTCQFTHAENGTNYIVTATRSFTVEEKGFSMFPVITNPPYNHTMEVEIGKPASIACSACFGKG SHFLADVLWQINKTVVGNFGEARIQEEEGRNESSNDMDCLTSVLRTGVTEKDLSLEYDCLALNLHGMIRHTI RLRRKQPSKECPSHIAVDHHHHHH
Formulation	Lyophilized from a 0.2 µm filtered solution of PBS, pH 7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7C for 2-7 days. Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only