

Recombinant Mouse IL-12 Receptor Subunit β2,IL-12RB2 (C-Fc)



Catalog No: #AP76208

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP76208-1 10ug #AP76208-2 50ug #AP76208-3 500ug #AP76208-4 1mg

Description

Product Name	Recombinant Mouse IL-12 Receptor Subunit β2,IL-12RB2 (C-Fc)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Mouse
Accession No.	P97378
Uniprot	P97378
GenID	16162;
Calculated MW	98.4kDa
Target Sequence	MAQTVRECSLALLFLFMWLLIKANIDVCKLGTVTVPAPVIPLGSAANISCSLNPKGCSHYPPSSNELILLKFVN DVLVENLHGKKVHDHTGHSSTFQVTNLSLGMTLFVCKLNCNSQKPPVPCVGVVEISVGVAPPEPPQNISCVQ EGENGTVACSWNSGKVTYLKTNYTLQLSGPNNLTCQKQCFSDNRQNCNRLDLGINLSPDLAESRFIVRVTA DLGNSSSLPHTFTFLDIVIPLPPWDIRINFLNASGSRGTLQWEDEGQVVLNQLRYQPLNSTSWNMVNATNAKG KYDLRDLRPFTEYEFQISSKLHLSGGWSNWSSESLRTRTPEEPPVGLDIWYMKQDIDYDRQQISLFWKSLNP SEARGKILHYQVTLQEVTKKTTLQNTTRHTSWTRVIPRTGAWTASVSAANSKGASAPTHINIVDLCGTGLLAPH QVSAKSENMDNILVTWQPPKADSAVREYIVWRALQPGSITKFPFHWLRIPPDNMSALISENIKPYICYEIRVH ALSESQGGCSSIRGDSKHKAPVSGPHITAITEKKERLFIWTHIPFPEQRGCILHYRIYWKERDSTAQPELCEIQ YRRSQNSHPISLQPRVTVLWMTAVTAAGESPQGNREFCPQGANVDDIEGRMDEPKSCDKTHTCPPCP APELLGGPSVFLFPPKPKDTLMISRTPEVTCVVVDVSHEDPEVKFNWYVDGVEVHNAKTKPREEQYNSTYRV VSVLTVLHQDWLNGKEYKCKVSNKALPAPIEKTISKAKGQPREPQVYTLPPSREEMTKNQVSLTCLVKGFYPS DIAVEWESNGQPENNYKTPPVLDSDGSFFLYSKLTVDKSRWQQGNVFCSSVMHEALHNHYTQKSLSLSPG K
Formulation	Lyophilized from a 0.2 μm filtered solution of PBS, pH 7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7C for 2-7 days. Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only