

Recombinant Human Peroxiredoxin-5, PRDX5 (N-6His)

Catalog No: #AP76209

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP76209-1 10ug #AP76209-2 50ug #AP76209-3 500ug #AP76209-4 1mg

Description

Product Name	Recombinant Human Peroxiredoxin-5, PRDX5 (N-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P30044
Uniprot	P30044
GeneID	25824;
Calculated MW	17.9kDa
Target Sequence	HHHHHHMAPIKVGDAIPAVEVFEGEPGNKVNLAELFKGKKGVLFGVPGAFTPGCSKTHLPGFVEQAEALKAK GVQVVAACLSVNDAFVTGEWGRAHKAEGKVRLADPTGAFGKETDLLLDDSLVSIFGNRRLLKRFSMVVQDGIV KALNVEPDGTGLTCSLAPNIISQL
Formulation	Lyophilized from a 0.2 µm filtered solution of PBS, pH 7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7C for 2-7 days. Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only