

## Recombinant Mouse CD83,HB15 (C-Fc)

Catalog No: #AP76214



Package Size: #AP76214-1 10ug #AP76214-2 50ug #AP76214-3 500ug #AP76214-4 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Mouse CD83,HB15 (C-Fc)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Mouse
Accession No.	O88324
Uniprot	O88324
GeneID	12522;
Calculated MW	39.5kDa
Target Sequence	MREVTVACSETADLPCTAPWDPQLSYAVSWAKVSESGTESVELPESKQNSSFEAPRRRRAYSLSLIQNTTICSS GTYRCALQELGGQRNLSGTVVLKVTGCPKEATESTFRKYRAVDDIEGRMDEPKSCDKTHTCPPCPAPPELLGG PSVFLFPPKPKDTLMISRTPEVTCVVVDVSHEDPEVKFNWYVDGVEVHNAKTKPREEQYNSTYRVVSVLTVLH QDWLNGKEYKCKVSNKALPAPIEKTISKAKGQPREPQVYTLPPSREEMTKNQVSLTCLVKGFYPSDIAVEWES NGQPENNYKTTTPVLDSDGSFFLYSKLTVDKSRWQQGNVFCFSVMHEALHNHYTQKSLSLSPGK
Formulation	Lyophilized from a 0.2 µm filtered solution of PBS,pH7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only