

Recombinant Mouse Growth Hormone Receptor, GHR, GHBP (C-Fc)



Catalog No: #AP76215

Orders: order@signalwayantibody.com

Package Size: #AP76215-1 10ug #AP76215-2 50ug #AP76215-3 500ug #AP76215-4 1mg

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Mouse Growth Hormone Receptor, GHR, GHBP (C-Fc)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Mouse
Accession No.	Q3UP14
Uniprot	Q3UP14
Calculated MW	58.3kDa
Target Sequence	MDLCQVFLTLALAVTSSTFSGSEATPATLGKASPVLQRINPSLGTSSSGKPRFTKCRSPELETFSCYWTEGDN PDLKTPGSIQLYYAKRESQRQAARIAHEWTQEWKECPDYVSAGKNSCYFNSSYTSIWIPYCIKLTNGDLLDQ KCFTVDEIVQPDPPIGLNWLLNISLTGIRGDIQVSWQPPNADV LKGWIILEYEIQYKEVNESKWKVMGPIWLT YCPVYSLRMDKEHEVRVRSRQRSFEKYSEFSEVLRVIFPQTNILEACEEDIQVDDIEGRMDEPKSCDKTHTCP PCPAPELLGGPSVFLFPPKPKDTLMISRTPEVTCVVVDVSHEDPEVKFNWYVDGVEVHNAKTKPREEQYNST YRVVSVLTVLHQDWLNGKEYKCKVSNKALPAPIEKTISKAKGQPREPQVYTLPPSREEMTKNQVSLTCLVKGF YPSDIAVEWESNGQPENNYKTTTPVLDSDGSFFLYSKLTVDKSRWQQGNV FSCVMHEALHNHYTQKSLSLS PGK
Formulation	Lyophilized from a 0.2 µm filtered solution of PBS, pH7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7C for 2-7 days. Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only