

Recombinant Mouse Transforming Growth Factor- β Receptor Type II, TGFBR2 (C-6His)

Catalog No: #AP76218

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP76218-1 10ug #AP76218-2 50ug #AP76218-3 500ug #AP76218-4 1mg

Description

Product Name	Recombinant Mouse Transforming Growth Factor- β Receptor Type II, TGFBR2 (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Mouse
Accession No.	Q62312-2
Uniprot	Q62312
GeneID	21813;
Calculated MW	16.2kDa
Target Sequence	IPPHVPKSVNSDVMASDNGGAVKLPQLCKFCDVRLSTCDNQKSCMSNCSITAICEKPHEVCVAVWRKNDKNIT LETVCHDPKLYHGFTLEDAASPKCVMKEKKRAGETFFMCACNMEECNDYIIFSEEYTTSSPDVDHHHHHH
Formulation	Lyophilized from a 0.2 μ m filtered solution of PBS, pH7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7C for 2-7 days. Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only