

Recombinant Mouse Serpin F2,Alpha-2-Antiplasmin (C-Fc)



Catalog No: #AP76219

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP76219-1 10ug #AP76219-2 50ug #AP76219-3 500ug #AP76219-4 1mg

Description

Product Name	Recombinant Mouse Serpin F2,Alpha-2-Antiplasmin (C-Fc)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Mouse
Accession No.	Q61247
Uniprot	Q61247
GeneID	18816;
Calculated MW	79.3kDa
Target Sequence	VDLPGQQPVSEQAQQKLPALFKLDNQDFGDHATLKRSPGHCKSVPTAEETRRLAQAMMAFTTDLFSLVAQ TSTSSNLVLSPLSVALALSHLALGAQNQTLHSLHRVLMNTGSCLPHLLSHFYQNLGPGTIRLAARIYLQKGFPI KDDFLEQSERLFGAKPVKLTGKQEEDLANINQWVKEATEGKIEDFLSELPDSTVLLLLNAIHFHGFWRTKFDPS LTQKDFHLDERFTVSDMMHAVSYPLRWFLLEQPEIQVAHFPPKNNMSFVVVMPTFYFEWNVSEVLANTLWD TLYHPSLQERPTKVWLPKLHLQQQLDLVATLSQLGLQELFQGPDLRGISEQNLVSSVQHQSTMELSEAGVE AAAATSVAMNRMSLSSFTVNRPFLFFIMEDTIGVPLFVGSVRNPNPSALPQLQEQRDSDPNRLIGQNDKADFH GGKTFGPDCLKAPRMEEDYPQFSSPKVDDIEGRMDEPKSCDKTHTCPPCPAPELLGGPSVFLFPPKPKDTLMI SRTPEVTCVVVDVSHEDPEVKFNWYVDGVEVHNAKTKPREEQYNSTYRVVSVLTVLHQDWLNGKEYKCKVVS NKALPAPIEKTKAKGQPREPQVYTLPPSREEMTKNQVSLTCLVKGFYPSDIAVEWESNGQPENNYKTTTPPV LDSDGSFFLYSKLTVDKSRWQQGNVFCFSVMHEALHNHYTQKSLSLSPGK
Formulation	Supplied as a 0.2 µm filtered solution of 20mM Tris,150mM NaCl,pH 8.0(20% glycerol,1mM DTT).
Storage	Store at < -20C, stable for 6 months after receipt.Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only