

# Recombinant Mouse Serpin F2,Alpha-2-Antiplasmin (C-6His)



Catalog No: #AP76220

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Package Size: #AP76220-1 10ug #AP76220-2 50ug #AP76220-3 500ug #AP76220-4 1mg

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Mouse Serpin F2,Alpha-2-Antiplasmin (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Mouse
Accession No.	Q61247
Uniprot	Q61247
GeneID	18816;
Calculated MW	53.2kDa
Target Sequence	VDLPGQQPVSEQAQQKLPALFKLDNQDFGDHATLKRSPGHCKSVPTAEETRRLAQAMMAFTTDLFSLVAQ TSTSSNLVLSPLSVALALSHLALGAQNQTLHSLHRVLMNTGSCPLHLLSHFYQNLGPGTIRLAARIYLQKGFPI KDDFLEQSERLFGAKPVKLTGKQEEDLANINQWVKEATEGKIEDFLSELPDSTVLLLLNAIHFHGFWRTKFDPS LTQKDFHLDERFTVSDMMHAVSYPLRWFLLEQPEIQVAHFPPKNNMSFVVVMPTYFEWNVSEVLANTWD TLYHPSLQERPTKVWPKLHLQQQLDLVATLSQLGLQELFQGPDLRGISEQNLVSSVQHQMELSEAGVE AAAATSVAMNRMSLSSFTVNRPFLLFIMEDTIGVPLFVGSVRNPNPSALPQLQEQRDSDPNRLIGQNDKADFH GGKTFGPDCLKLAPRMEEDYPQFSSPKVDHHHHHH
Formulation	Supplied as a 0.2 µm filtered solution of 20mM Tris, 150mM NaCl, pH 8.0 (20% glycerol, 1mM DTT).
Storage	Store at < -20C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only