

Recombinant Mouse Tissue Inhibitors of Metalloproteinases 1,TIMP1

Catalog No: #AP76223

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP76223-1 10ug #AP76223-2 50ug #AP76223-3 500ug #AP76223-4 1mg

Description

Product Name	Recombinant Mouse Tissue Inhibitors of Metalloproteinases 1,TIMP1
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Mouse
Accession No.	P12032
Uniprot	P12032
GeneID	21857;
Calculated MW	20.2kDa
Target Sequence	CSCAPPHPQTAF CNSDLVIRAKFMGSPEINETTLYQRYKIKMTKMLKGFKA VGNAA DIRYAYTPVMESLCGYA HKSQNRSEEF LITGRLRNGNLHISACSF LVPWRTLSPAQQRAFSKTY SAGCGVCTVFPCL SIPCKLES DTHCL WTDQVLVGS EDYQSRHFACLPRNPGLCTWRSLGAR
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM Tris, 150mM NaCl, pH8.0.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7C for 2-7 days. Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only