

Recombinant Mouse NGAL,Lipocalin-2,LCN2 (C-6His)

Catalog No: #AP76229



Package Size: #AP76229-1 10ug #AP76229-2 50ug #AP76229-3 500ug #AP76229-4 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Mouse NGAL,Lipocalin-2,LCN2 (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Mouse
Accession No.	P11672
Uniprot	P11672
GeneID	16819;
Calculated MW	21.9kDa
Target Sequence	QDSTQNLIPAPSLTVPPLQPDFRSDQFRGRWYVGLAGNAVQKKTEGSFTMYSTIYELQENNSYNVTSILVRD QDQGCRYWIRTFVPSSRAGQFTLGNMHRYPQVQSYNVQVATTDYNQFAMVFFRKTSENKQYFKITLYGRTK ELSPELKERFTRFAKSLGLKDDNIIFSVPTDQCIDNVDDHHHHH
Formulation	Supplied as a 0.2 µm filtered solution of 20mM MES, 150mM NaCl, 10% Glycerol, pH 5.5.
Storage	Store at < -20C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only