

Recombinant Mouse NGAL,Lipocalin-2,LCN2 (C-Fc)

Catalog No: #AP76230



Package Size: #AP76230-1 10ug #AP76230-2 50ug #AP76230-3 500ug #AP76230-4 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Mouse NGAL,Lipocalin-2,LCN2 (C-Fc)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Mouse
Accession No.	P11672
Uniprot	P11672
GeneID	16819;
Calculated MW	48kDa
Target Sequence	<p>QDSTQNLIPAPSLLTVP LQPDRSDQFRGRWYVGLAGNAVQKKTEGSFTMYSTIYELQENNSYNVTSILVRD</p> <p>QDQGCRYWIRTFVPSSRAGQFTLGNMHRYPQVQSYNVQVATTDYDYNQFAMVFFRKTSENKQYFKITLYGRTK</p> <p>ELSPELKERFTRFAKSLGLKDDNIIFSVPTDQCIDNVDDIEGRMDEPKSCDKTHTCPPCPAPPELLGGPSVFLFP</p> <p>PKPKDTLMISRTPEVTCVVVDVSHEDPEVKFNWYVDGVEVHNAKTKPREEQYNSTYRVVSVLTVLHQDWLNG</p> <p>KEYKCKVSNKALPAPIEKTISKAKGQPREPQVYTLPPSREEMTKNQVSLTCLVKGFYPSDIAVEWESNGQPEN</p> <p>NYKTTTPVLDSDGSFFLYSKLTVDKSRWQQGNVFCFSVMHEALHNHYTQKSLSLSPGK</p>
Formulation	Supplied as a 0.2 µm filtered solution of 20mM MES, 150mM NaCl, 10% Glycerol, pH 5.5.
Storage	Store at < -20C, stable for 6 months after receipt.Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only