

# Recombinant Rat Fibroblast Growth Factor 2,FGF2,FGFb



Catalog No: #AP76233

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Package Size: #AP76233-1 10ug #AP76233-2 50ug #AP76233-3 500ug #AP76233-4 1mg

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Rat Fibroblast Growth Factor 2,FGF2,FGFb
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Rat
Accession No.	P13109
Uniprot	P13109
GeneID	54250;
Calculated MW	16.2kDa
Target Sequence	ALPEDGGGAFPPGHFKDPKRLYCKNGGFFLRHPDGRVDGVREKSDPHVKLQLQAEERGVSISIKGVCANRYL AMKEDGRLLASKCVTEECFFFERLESNNYNTYRSRKYSWYVALKRTGQYKLGSKTGPGQKAILFLPMSAKS
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM PB, 150mM NaCl,pH7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only