

Recombinant Human Phosphoglucomutase 2,PGM2 (N-6His)



Catalog No: #AP76234

Orders: order@signalwayantibody.com

Package Size: #AP76234-1 10ug #AP76234-2 50ug #AP76234-3 500ug #AP76234-4 1mg

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Phosphoglucomutase 2,PGM2 (N-6His)
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q96G03
Uniprot	Q96G03
GeneID	55276;
Calculated MW	70.5kDa
Target Sequence	MGSSHHHHHHSSGLVPRGSHMAAPEGSGLDEDARLDQETAQWLRWDKNSLTLEAVKRLIAEGNKEELRKCFC GARMEFGTAGLRAAMGPGISRMNDLTIQTTQGFCRYLEKQFSDLKQKGVISFDARAHPSGGSSRRFARLA ATTFISQGIPVYLFSDITPTFPVFTVSHLKLKAGIMITASHNPKQDNGYKVYWDNGAQIISPHDKGISQAIENL EPWPQAWDDSLIDSSPLLHNPSASINNDYFEDLKKYCFHRSVNRETKVKFVHTSVHGVGHSFVQSAFKAFDL VPPEAVPEQKDPDPEFPTVKYPNPEEGKGVLTLSFALADKTKARIVLANDPADRLAVAEKQDSEGEWRVFSG NELGALLGWWLFTSWKEKNQDRSALKDYMVLSSTVSSKILRAIALKEGFHFEETLTGFKWWMGNRAKQLIDQG KTVLFAFEEAIGYMCCPFVLDKDGVSAAVISAEASFLATKNLSLQQLKAIYVEYGYHITKASYFICHQDQETIKK LFENLRNYDGKNNYPKACGKFEISAIRDLTTGYDDSQPKKAVLPTSKSSQMITFTFANGGVATMRTSGTEPKI KYAELCAPPGNSDPEQLKELNELVSAIEEHFFQPQKYNLQPKAD
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM Tris200mM NaClpH8.0.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only