

Recombinant Human Follicle-Stimulating Hormone α,β Dimer



Catalog No: #AP76235

Orders: order@signalwayantibody.com

Package Size: #AP76235-1 10ug #AP76235-2 50ug #AP76235-3 500ug #AP76235-4 1mg

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Follicle-Stimulating Hormone α,β Dimer
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P01215&P01225
Uniprot	P01215
GeneID	1081;
Calculated MW	10.2&12.5kDa
Target Sequence	APDVQDCPECTLQENPFFSQPGAPILQCMGCCFSRAYPTPLRSKKTMLVQKNVTSESTCCVAKSYNRVTVM GGFKVENHTACHCSTCYHKS&NSCELNTNIAIEKEECRFCISINTTWCAGYCYTRDLVYKDPARPKIQKTCTF KELVYETVRVPGCAHHADSLYTPVATQCHCGKCDSSTDCTVRGLGPSYCSFGEMKE
Formulation	Lyophilized from a 0.2 μ m filtered solution of PBS, pH7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only