

Recombinant Mouse Cathepsin S,CTSS (C-6His)

Catalog No: #AP76242



Package Size: #AP76242-1 10ug #AP76242-2 50ug #AP76242-3 500ug #AP76242-4 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Mouse Cathepsin S,CTSS (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Mouse
Accession No.	O70370
Uniprot	O70370
GeneID	13040;
Calculated MW	37.5kDa
Target Sequence	VCSVAMEQLQRDPTLDYHWDLWKKTHEKEYKDKNEEEVRRLLIWEKNLKFIMIHNLEYSMGMHYQVGMNDM GDMTNEEILCRMALRIPRQSPKTVTFRSYSNRTLPTVDWREKGCVTEVKYQGSCGACWAFSAVGALEGQ LKLKTGKLISLSAQNLVDCSNEEKYGNKGGGYMTEAFQYIIDNGGIEADASYPYKATDEKCHYNSKNRAAT CSRYIQLPFGDEDALKEAVATKGPVSVGIDASHSSFFFYKSGVYDDPSCTGNVNHGVLVVGYGTLDGKDYWL VKNSWGLNFGDQGYIRMARNNKNHCGIASYCSYPEIVDHHHHHH
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM PB,150mM NaCl,pH7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only