

Recombinant Mouse Interleukin-12 Subunit β ,IL-12 p40,IL-12B (C-6His)



Catalog No: #AP76243

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP76243-1 10ug #AP76243-2 50ug #AP76243-3 500ug #AP76243-4 1mg

Description

Product Name	Recombinant Mouse Interleukin-12 Subunit β ,IL-12 p40,IL-12B (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Mouse
Accession No.	P43432
Uniprot	P43432
GeneID	16160;
Calculated MW	37.1kDa
Target Sequence	MWELEKDVYVVEVDWTPDAPGETVNLTCDTPEEDDITWTSQQRHGVIGSGKTLTITVKEFLDAGQYTCHKGG ETLSHSHLLLHKKENGIWSTEILKNFKNKTKFLKCEAPNYSGRFTCSWLVQRNMDLKFNIKSSSSSPDSRAVTC GMASLSAEKVTLQDRDYEKYSVSCQEDVTCPTAEETLPIELALEARQQNKYENYSTSFFIRDIKPDPPKNLQM KPLKNSQVEVSWEYPDSWSTPHSYFSLKFFVRIQRKKEKMKETEEGCNQKGAFLVEKTSTEVQCKGGNVCV QAQDRYYNSSCSKWACVPCRVRSDHHHHHHH
Formulation	Lyophilized from a 0.2 μ m filtered solution of PBS, pH7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only