

Recombinant Mouse Interleukin-12,IL-12

Catalog No: #AP76244



Package Size: #AP76244-1 10ug #AP76244-2 50ug #AP76244-3 500ug #AP76244-4 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Mouse Interleukin-12,IL-12
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Mouse
Accession No.	P43432 & NP032377
Uniprot	P43432
GeneID	16160;
Calculated MW	57.5kDa
Target Sequence	MWELEKDVYVVEVDWTPDAPGETVNLTCDTPEEDDITWTSQQRHGVIGSGKTLTITVKEFLDAGQYTCHKGG ETLSHSHLLLHKKENGIWSTEILKNFKNKTKFLKCEAPNYSGRFTCSWLVRNMDLKFNIKSSSSSPDSRAVTC GMASLSAEKVTLQQRDYEKYSVSCQEDVTCPTAEETLPIELALEARQQNKYENYSTSFFIRDIKPDPPKLNLM KPLKNSQVEVSWEYPDSWSTPHSYFSLKFFVRIQRKKEKMKETEEGCNQKGAFLVEKTSTEVQCKGGNVCV QAQDRYNNSSCSKWACVPCRVRS&RVIPVSGPARCLSQSRNLLKTTDDMVKTAREKLKHYSCTAEDIDHEDIT RDQTSTLKTCLPELHKNESCLATRETSSTTRGSCLPQKTSMMMLCLGSIYEDLKMVQTEFQAINAALQNH HQQIILDKGMLVAIDELMQSLNHNGETLRQKPPVGEADPYRVKMKLCILLHAFSTRVVTINRVMGYLSSA
Formulation	Lyophilized from a 0.2 µm filtered solution of PBS, pH7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only