

# Recombinant Mouse Interferon $\zeta$ , IFN- $\zeta$ , Limitin (C-6His)

Catalog No: #AP76250

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Package Size: #AP76250-1 10ug #AP76250-2 50ug #AP76250-3 500ug #AP76250-4 1mg

## Description

Product Name	Recombinant Mouse Interferon $\zeta$ , IFN- $\zeta$ , Limitin (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Mouse
Accession No.	Q9R1T0
Uniprot	Q9R1T0
Calculated MW	19.5kDa
Target Sequence	LDSGKSGSLHLERSETARFLAELRSVPGHQCLRDRTDFPCPWKEGTNITQMTLGETTSCYSQTLRQVLHLFDT EASRAAWHERALDQLLSSLWRELQVLKSPREQQGQSCPLPFALAIRTYFRGFFRYLKAKAHSACSWEIVRVQL QVDLPAFPLSARRGPRVDHHHHHH
Formulation	Lyophilized from a 0.2 $\mu$ m filtered solution of PBS, pH7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7C for 2-7 days. Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only