

## Recombinant Mouse DDR2,CD167,TKT (C-6His)

Catalog No: #AP76252



Package Size: #AP76252-1 10ug #AP76252-2 50ug #AP76252-3 500ug #AP76252-4 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Mouse DDR2,CD167,TKT (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Mouse
Accession No.	Q62371
Uniprot	Q62371
GeneID	18214;
Calculated MW	43.5kDa
Target Sequence	QVNPAICRYPLGMSGGHIPDEITASSQWSESTAAKYGRLDSEEGDGAWCPEIPVQPDDLKEFLQIDLRTLHFI TLVGTQGRHAGGHGIEFAPMYKINYSRDGSRWISWRNRHGKQVLDGNSNPYDVFLKDLEPPIVARFVRLIPVT DHSMNVCMRVELYGCVWLDGLVSYNAPAGQQFVLPGGSIILNDSVYDGAVGYSMTEGLGQLTDGVSGLDD FTQTHEYHVWPGYDYVGWRNESATNGFIEIMFEFDRIRNFTTMKVHCNMFAGVKIFKEVQCYFRSEASEW EPTAVYFPLVLDVNP SARFVTVPLHHRMASAIKCQYHFADTWMMFSEITFQSDAAMYNNSGALPTSPMAPT TYDPMLKVDDSNTRVDHHHHHHH
Formulation	Lyophilized from a 0.2 µm filtered solution of PBS, pH7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only