

Recombinant Human Mesothelin,MSLN,CAK1,MPF (C-Fc)



Catalog No: #AP76255

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP76255-1 10ug #AP76255-2 50ug #AP76255-3 500ug #AP76255-4 1mg

Description

Product Name	Recombinant Human Mesothelin,MSLN,CAK1,MPF (C-Fc)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q13421
Uniprot	Q13421
GeneID	10232;
Calculated MW	53.9kDa
Target Sequence	LAGETGQEAAPLDGVLNPPNISSLSPRQLLGFPCEVSGLSTERVRELAVALAQKNVKLSTEQLRCLAHRLS EPPEDLDALPLDLLLFLNPDAFSGPQACTRFFSRITKANVDLLPRGAPERQRLLPAALACWGVRSLLSEADV RALGGLACDLPGRFVAESAEVLLPRLVSCPGPLDQDQQAARAALQGGGPPYGGPSTWSVSTMDALRGLLP VLGQPIIRSIPQGIVAAWRQRSSRDPSWRQPERVDDIEGRMDEPKSCDKTHTCPPCPAPELLGGPSVFLFPPK PKDTLMISRTPEVTCVVVDVSHEDPEVKFNWYVDGVEVHNAKTKPREEQYNSTYRVVSVLTVLHQDWLNGKE YKCKVSNKALPAPIEKTISKAKGQPREPQVYTLPPSREEMTKNQVSLTCLVKGFYPSDIAVEWESNGQPENNY KTTTPVLDSDGSFFLYSKLTVDKSRWQQGNVFSCSVMHEALHNHYTQKSLSLSPGK
Formulation	Lyophilized from a 0.2 µm filtered solution of PBS, pH7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only