

Recombinant Mouse Follistatin-Like 1, FSTL1 (C-Fc)

Catalog No: #AP76260



Package Size: #AP76260-1 10ug #AP76260-2 50ug #AP76260-3 500ug #AP76260-4 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Mouse Follistatin-Like 1, FSTL1 (C-Fc)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Mouse
Accession No.	Q62356
Uniprot	Q62356
GeneID	14314;
Calculated MW	59.6kDa
Target Sequence	<p>EEEPRSKKICANVFCGAGRECAVTEKGEPTCLCIEQCKPHKRPVCGSNGKTYLNHCELHRDACTGSKIQVD</p> <p>YDGHCKEKKASPSASPVVICYQANRDELRRRLIQWLEAEIIPDGWFSKGSNYSEILDKYFKSFDNGDSHLDSS</p> <p>EFLKFVEQNETAINITTYADQENKLLRSLCVDALIELSDENADWKLFSQEFLKCLNPSFNPEKKCALEDETYA</p> <p>DGAETEVDNRCVCSCGHVWCTAMTCDGKNQKGVQTHTEEEKTGYVQELQKHGTAEKTKKVNTKEIVDDI</p> <p>EGRMDEPKSCDKTHTCPPCPAPPELLGGPSVFLFPPKPKDTLMISRTPEVTCVVVDVSHEDPEVKFNWYVDGV</p> <p>EVHNAKTKPREEQYNSTYRVVSVLTVLHQDWLNGKEYKCKVSNKALPAPIEKTISKAKGQPREPQVYTLPPSR</p> <p>EEMTKNQVSLTCLVKGFYPSDIAVEWESNGQPENNYKTTTPVLDSDGSFFLYSKLTVDKSRWQQGNVFCSSV</p> <p>MHEALHNHYTQKSLSLSPGK</p>
Formulation	Lyophilized from a 0.2 µm filtered solution of PBS, pH7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7C for 2-7 days. Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only