

Recombinant Mouse Prolactin Receptor,PRLR (C-Fc)

Catalog No: #AP76261



Package Size: #AP76261-1 10ug #AP76261-2 50ug #AP76261-3 500ug #AP76261-4 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Mouse Prolactin Receptor,PRLR (C-Fc)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Mouse
Accession No.	Q08501
Uniprot	Q08501
GeneID	19116;
Calculated MW	51.7kDa
Target Sequence	<p>QSPPGKPEIHKCRSPDKETFTCWNNPGSDGGLPTNYSLTYSKEGEKNTYECPDYKTS GPNSCFFSKQYTSI</p> <p>WKIYIITVNATNEMGSSTSDPLYVDVTYIVEPEPPRNLTEVKQLKDKKTYLWVKWLPPTITDVKTGWFTMEYEI</p> <p>RLKSEEADEWEIHFTGHQTQFKVFDLYPGQKYL VQTRCKPDHGYWSRWGQEKSIEIPNDFTLKDVDDIEGRM</p> <p>DEPKSCDKTHTCPPCPAPELLGGPSVFLFPPKPKDTLMISRTPEVTCVVDVSHEDPEVKFNWYVDGVEVHN</p> <p>AKTKPREEQYNSTYRVVSVLTVLHQDWLNGKEYKCKVSNKALPAPIEKTISKAKGQPREPQVYTLPPSREEMT</p> <p>KNQVSLTCLVKGFYPSDIAVEWESNGQPENNYKTPPVLDSDGSFFLYSKLTVDKSRWQQGNV FSCSVMHE</p> <p>ALHNHYTQKSLSLSPGK</p>
Formulation	Lyophilized from a 0.2 µm filtered solution of PBS, pH7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only