

Recombinant Mouse IAP,OA3,CD47 (C-Fc)

Catalog No: #AP76265



Package Size: #AP76265-1 10ug #AP76265-2 50ug #AP76265-3 500ug #AP76265-4 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Mouse IAP,OA3,CD47 (C-Fc)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Mouse
Accession No.	Q61735-2
Uniprot	Q61735
GeneID	16423;
Calculated MW	42.8kDa
Target Sequence	<p>QLLFSNVNSIEFTSCNETWIPCIVRNVEAQSTEEMFVKWKLNKSIFIFIYDGNKNSTTTDQNFSAKISVSDLING IASLKMDKRDAMVGNITCEVTELSREGKTVIELKNRTAFNTDQGSACSYYYYEKGCKLVSWFSPVDDIEGRM DEPKSCDKTHTCPPCPAPPELLGGPSVFLFPPKPKDTLMISRTPEVTCVVDVSHEDPEVKFNWYVDGVEVHN AKTKPREEQYNSTYRVVSVLTVLHQDWLNGKEYKCKVSNKALPAPIEKTISKAKGQPREPQVYTLPPSREEMT KNQVSLTCLVKGFPYPSDIAVEWESNGQPENNYKTTTPVLDSDGSFFLYSKLTVDKSRWQQGNVFSCSVMHE ALHNHYTQKSLSLSPGK</p>
Formulation	Lyophilized from a 0.2 µm filtered solution of PBS, pH7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7C for 2-7 days. Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only